

JointProject

Department of Cardiovascular Surgery  
Graduate School of Medicine  
Osaka University


×  
SANPLATEC®



A successful series of  
“live transport and display”  
events of myocardial sheets

with iP-TEC®

iPS cell-derived  
myocardial sheet



Major successful events for the transport  
and display of myocardial sheets

2016

May: Ise-Shima Summit  
Osaka ↔ Mie (Car)

September: Hyogo-Kobe Medical and Health Care  
Fair, G7 Kobe Health Ministers'  
Meeting

Osaka ↔ Kobe (Car)

October: Regenerative Medicine Japan  
Osaka ↔ Yokohama (Japan Railways)

2017

June: International Society for Stem Cell  
Research (ISSCR)  
Osaka ↔ San Francisco ↔ Boston  
(Airplane)

November: Stem Cell Society Singapore (SCSS)  
Osaka ↔ Malaysia ↔ Singapore  
(Airplane)

Watch how the myocardial  
sheet pulsates by scanning  
the QR code on the right.



iP-TEC® successfully offsets the harsh  
vibrations associated with baggage handling



While returning from overseas exhibitions (Boston, Singapore), we analyzed the vibration data from the logger attached to the tertiary container to find that the container had gone through harsh impacts and vibrations during loading and unloading from the airplane. However, thanks to iP-TEC®'s technology, the primary container is filled with culture solution without any bubbles inside, and the iP-TEC® secondary container stabilizes and retains the primary container, leaving the myocardial sheet intact.

In addition, we analyzed the temperature logger attached to inside of the secondary container and found out that the temperature-controlled transport box, composed of iP-TEC® heat storage materials and an iP-TEC® Premier BOX V-8.5, had maintained temperature very stably and contributed to successful transportation.